

Pricing is based on active, level of effort, and compound form and is dependent on Service Level Agreements; the tier placement below is for Injectable Solutions only.

Active	Method Type	Stability Indicating**	Standard Type	Method Number	Target Stability***	Minimum Amount Required****	Tier
5-Fluorouracil	HPLC/ UHPLC	Yes	USP	4017	30	5mL	Tier 2
7-keto DHEA	HPLC/ UHPLC	Yes	Medisca	4258	180	5mL	Tier 1
Acetaminophen	HPLC/ UHPLC	Yes	USP	4406	N/A	5mL	Tier 2
Adenosine	HPLC/ UHPLC	Yes	USP	4319	45	5mL	Tier 1
Alprostadil	HPLC/ UHPLC	Yes	USP	4232	180	5mL	Tier 2
Amiodarone HCl	HPLC/ UHPLC	Yes	USP	4335	120	5mL	Tier 1
Amiodarone	HPLC/ UHPLC	Yes	USP	4335	120	5mL	Tier 1
Amitriptyline HCl	HPLC/ UHPLC	Yes	USP	4269	90	5mL	Tier 1
Amitriptyline	HPLC/ UHPLC	Yes	USP	4269	91	5mL	Tier 1
Amlodipine Besylate	HPLC/ UHPLC	Yes	USP	4357	N/A	5mL	Tier 2
Ampicillin	HPLC/ UHPLC	Yes	USP	4267	160	5mL	Tier 3
Anastrozole	HPLC/ UHPLC	Yes	USP	4264	180	5mL	Tier 1
Apomorphine	HPLC/ UHPLC	Yes	USP	4448	N/A	5mL	Tier 2
Ascorbic Acid	HPLC/ UHPLC	Yes	USP	4350	180	5mL	Tier 2
Atropine Sulfate	HPLC/ UHPLC	Yes	USP	4304	90	5mL	Tier 1
Atropine	HPLC/ UHPLC	Yes	USP	4304	90	5mL	Tier 1
Baclofen	HPLC/ UHPLC	Yes	USP	4010	180	5mL	Tier 1
Benzocaine	HPLC/ UHPLC	Yes	USP	4328	60	5mL	Tier 1
Benzyl Alcohol	HPLC/ UHPLC	Yes	USP	4409	N/A	5mL	Tier 2
Betamethasone (Base)	HPLC/ UHPLC	Yes	USP	4278	90	5mL	Tier 1
Betamethasone Acetate	HPLC/ UHPLC	Yes	USP	4202	90	5mL	Tier 2
Betamethasone Sodium Phosphate	HPLC/ UHPLC	Yes	USP	4203	90	5mL	Tier 1
Biotin	HPLC/ UHPLC	Yes	USP	4341	180	5mL	Tier 1
Budesonide	HPLC/ UHPLC	Yes	USP	4385	N/A	5mL	Tier 2
Bumetanide	HPLC/ UHPLC	Yes	USP	4201	90	5mL	Tier 1
Bupivacaine HCl	HPLC/ UHPLC	Yes	USP	4399	90	5mL	Tier 1
Bupivacaine	HPLC/ UHPLC	Yes	USP	4399	90	5mL	Tier 1
Buprenorphine HCl	HPLC/ UHPLC	Yes	USP	4072	90	5mL	Tier 2

Pricing is based on active, level of effort, and compound form and is dependent on Service Level Agreements; the tier placement below is for Injectable Solutions only.

Active	Method Type	Stability Indicating**	Standard Type	Method Number	Target Stability***	Minimum Amount Required****	Tier
Buprenorphine	HPLC/ UHPLC	Yes	USP	4072	90	5mL	Tier 2
Caffeine Citrate	HPLC/ UHPLC	Yes	USP	4330	60	5mL	Tier 1
Caffeine	HPLC/ UHPLC	Yes	USP	4330	75	5mL	Tier 1
Calcium Chloride Dihydrate	IC	No	USP	6001	120	5mL	Tier 3
Calcium Citrate	IC	No	USP	6001	120	5mL	Tier 2
Calcium Gluconate	Titration	No	Volumetric Solution	5018	120	Nominally Equivalent to 2 gram of Calcium Gluconate	Tier 3
Carbamazepine	HPLC/ UHPLC	Yes	USP	4299	N/A	5mL	Tier 1
Cefazolin Sodium	HPLC/ UHPLC	Yes	USP	4331	90	5mL	Tier 1
Cefazolin	HPLC/ UHPLC	Yes	USP	4331	90	5mL	Tier 1
Cefepime	HPLC/ UHPLC	Yes	USP	4333	90	5mL	Tier 1
Ceftazidime	HPLC/ UHPLC	Yes	USP	4256	180	5mL	Tier 2
Ceftriaxone Sodium	HPLC/ UHPLC	Yes	USP	4101	100	5mL	Tier 3
Ceftriaxone	HPLC/ UHPLC	Yes	USP	4101	100	5mL	Tier 3
Cefuroxime	HPLC/ UHPLC	Yes	USP	4300	180	5mL	Tier 1
Cefuroxime Sodium	HPLC/ UHPLC	Yes	USP	4300	180	5mL	Tier 1
Choline Chloride	IC	Needs Quote	Medisca	6004	180	5mL	
Ciprofloxacin HCl	HPLC/ UHPLC	Yes	USP	4334	N/A	5mL	Tier 1
Ciprofloxacin	HPLC/ UHPLC	Yes	USP	4334	N/A	5mL	Tier 1
Clindamycin Phosphate	HPLC/ UHPLC	Yes	USP	4307	180	5mL	Tier 1
Clonidine HCl	HPLC/ UHPLC	Yes	USP	4253	90	5mL	Tier 1
Clonidine	HPLC/ UHPLC	Yes	USP	4253	90	5mL	Tier 1
Cyanocobalamin	HPLC/ UHPLC	Yes	Medisca	4040	180	5mL	Tier 1
Cyclobenzaprine HCl	HPLC/ UHPLC	Yes	USP	4249	120	5mL	Tier 1
Cyclobenzaprine	HPLC/ UHPLC	Yes	USP	4249	120	5mL	Tier 1

Pricing is based on active, level of effort, and compound form and is dependent on Service Level Agreements; the tier placement below is for Injectable Solutions only.

Active	Method Type	Stability Indicating**	Standard Type	Method Number	Target Stability***	Minimum Amount Required****	Tier
Cyclopentolate HCl	HPLC/ UHPLC	Yes	USP	4288	150	5mL	Tier 1
Cyclosporine	HPLC/ UHPLC	Yes	USP	4449	N/A	5mL	Tier 2
Dehydroepianandrosterone (DHEA)	HPLC/ UHPLC	Yes	Medisca	4233	180	5mL	Tier 1
Dexamethasone Acetate	HPLC/ UHPLC	Yes	USP	4055	90	5mL	Tier 2
Dexamethasone Sodium Phosphate	HPLC/ UHPLC	Yes	USP	4069	90	5mL	Tier 1
Dexamethasone Phosphate	HPLC/ UHPLC	Yes	USP	4069	90	5mL	Tier 1
Dexmedetomidine (Precedex)	HPLC/ UHPLC	Yes	USP	4052	90	5mL	Tier 2
Dextromethorphan Hydrobromide	HPLC/ UHPLC	Yes	USP	4282	150	5mL	Tier 2
Diazepam	HPLC/ UHPLC	Yes	USP	4247	90	5mL	Tier 1
Diclofenac	HPLC/ UHPLC	Yes	USP	4214	45	5mL	Tier 1
Diclofenac Potassium	HPLC/ UHPLC	Yes	USP	4214	N/A	5mL	Tier 1
Diclofenac Sodium	HPLC/ UHPLC	Yes	USP	4214	45	5mL	Tier 1
Diltiazem HCl	HPLC/ UHPLC	Yes	USP	4359	90	5mL	Tier 1
Diltiazem	HPLC/ UHPLC	Yes	USP	4359	90	5mL	Tier 1
Diphenhydramine HCl	HPLC/ UHPLC	Yes	USP	4376	50	5mL	Tier 2
Dipyridamole	HPLC/ UHPLC	Yes	USP	4291	60	5mL	Tier 1
Dobutamine	HPLC/ UHPLC	Yes	USP	4270	120	5mL	Tier 1
Dopamine HCl	HPLC/ UHPLC	Yes	USP	4393	90	5mL	Tier 1
Dopamine	HPLC/ UHPLC	Yes	USP	4393	90	5mL	Tier 1
Doxepin HCl	HPLC/ UHPLC	Yes	USP	4283	N/A	5mL	Tier 2
Doxepin	HPLC/ UHPLC	Yes	USP	4283	N/A	5mL	Tier 2
Doxycycline	HPLC/ UHPLC	Yes	USP	4420	N/A	5mL	Tier 2
Doxycycline Hyclate	HPLC/ UHPLC	Yes	USP	4420	N/A	5mL	Tier 2
Droperidol	HPLC/ UHPLC	Yes	USP	4329	75	5mL	Tier 1

Pricing is based on active, level of effort, and compound form and is dependent on Service Level Agreements; the tier placement below is for Injectable Solutions only.

Active	Method Type	Stability Indicating**	Standard Type	Method Number	Target Stability***	Minimum Amount Required****	Tier
Edetate Calcium Disodium	Titration	No	Volumetric Solution	5020	90	Nominally Equivalent to 2 gram of Edetate Calcium Disodium	Tier 2
Edetate Disodium	Titration	No	Volumetric Solution	5016 and 5017	180	Nominally Equivalent to 5 grams of Edetate Disodium	Tier 2
Ephedrine Sulfate	HPLC/ UHPLC	Yes	USP	4374	180	5mL	Tier 1
Ephedrine	HPLC/ UHPLC	Yes	USP	4374	180	5mL	Tier 1
Epinephrine Bitartrate	HPLC/ UHPLC	Yes	USP	4384	180	5mL	Tier 1
Epinephrine HCl	HPLC/ UHPLC	Yes	USP	4384	180	5mL	Tier 1
Epinephrine	HPLC/ UHPLC	Yes	USP	4384	180	5mL	Tier 1
Estradiol	HPLC/ UHPLC	Yes	USP	4209	540	5mL	Tier 1
Estradiol Valerate	HPLC/ UHPLC	Yes	USP	4349	180	5mL	Tier 2

Pricing is based on active, level of effort, and compound form and is dependent on Service Level Agreements; the tier placement below is for Injectable Solutions only.

Active	Method Type	Stability Indicating**	Standard Type	Method Number	Target Stability***	Minimum Amount Required****	Tier
Estriol	HPLC/ UHPLC	Yes	USP	4212	180	5mL	Tier 1
Estrone	HPLC/ UHPLC	Yes	USP	4221	180	5mL	Tier 1

*Method Feasibility has been confirmed and the Method Development is in progress, the target for completion of these methods is 2018. If you have a need for any of these actives sooner or any actives not listed, please email MarketingTechServices@dynalabs.us for Feasibility Quote or a Method Development Quote.

**Stability Indicating Methods are validated at the API level according to USP 1225 Category 1, complete with degradation studies. Validation Extension is still required for each formulation to meet the specificity requirement. For methodology that is not stability indicating due to limitations to separate degradants (e.g. IC and Titration) please contact MarketingTechServices@dynalabs.us for information about Stability Profile testing.

*****and Italics** Target Stability is not confirmed. Please email MarketingTechServices@dynalabs.us for more information and for suggested storage conditions.

****For all titrations the amount of sample required is concentration based and the amount of material listed is based on the amount of active required in some cases. Please email MarketingTechServices@dynalabs.us

This list is for information only; if there are any questions or concerns about the relevance of the data conveyed or the impact the information to your product please email MarketingTechServices@dynalabs.us with questions.